

Se bilder i bildarkiv

1 Dokumentbeskrivning

Dokumentet beskriver hur man kan se bilder i bildarkivet

1.1 Allmänt


Bilderna nås via ett uthopp i Cosmic *Bildarkiv röntgen/ klinfys* och via en länk i svaret på undersökningen.

2 Instruktioner

1. Logga in i Cosmic, välj aktuell patient.
2. Välj först Meny därefter *Beställning och Svar*
3. Välj därefter ”uthoppet”:

Bildarkiv Röntgen/ Klin fys

Följande fönster visas med bilder eller svar ifrån klinisk fysiologi.



The screenshot displays the SECTRA UniView software interface. At the top, the patient information bar shows: "Person, Namnlös, MR", "Född 1910-10-10", "Ålder 111 år", "Kön Man", and "Personnummer 1910101010". The user "Hedlund Jessami [jehed1] Administrator" is logged in. On the left, a sidebar contains icons for "Dokument", "Bilder", and "Ladda upp". The main area features a search bar "Filtrera undersökningar" and a list of examinations. The first examination is "EKO" dated "2022-09-02". Other examinations include "Nygranskn+Dem av Klin.Fys us" (2022-08-18), "VILO-EKG, datortolkning" (2022-08-16), and several "EKO" exams from August and July 2022. Some exams are expanded, showing "Inga bilder". The right side of the interface shows a Siemens CLV ANESTESIKLINIKEN header, technical parameters (VF13-5, C-Vascular, 30 dB, 8.9 MHz, DR 60 dB, Edge 2, Persist 2, R/S 3, Map F, Tint 1, 47 fps), and a large dark area for the image. A status bar at the bottom indicates "P 100% MI 0.5".