



Инструкция по сбору образца кала



Соберите кал.

Используйте бумагу для сбора или тщательно очищенный контейнер.




Отвинтите верхнюю часть маленькой пробирки для образца.

Под верхней частью прикреплена ложка.



Используя ложку, перенесите кал в небольшую пробирку для образца.


Заполните пробирку образцом кала на четверть.



Осторожно завинтите верхнюю часть маленькой пробирки для образца.



Напишите на **обеих** пробирках свой полный личный идентификационный номер и имя, а также дату и время сбора образца. Используйте этикетки, которые предоставил ваш доктор.



Поместите маленькую пробирку для образца в большую пробирку и закрутите её крышкой.



Как можно скорее доставьте образец своему врачу.

Сохраняйте пробирку с образцом охлаждённой всё время.

Имейте при себе действительное удостоверение личности.

Ответственное подразделение:

.....

По вопросам о сборе образца обращайтесь к:

.....

Адрес доставки образца:

.....

Clinical Microbiology for Kronoberg and Blekinge County

This information has been produced and prepared by medical laboratory scientist Ban Faraj and microbiologist Anna Pregiel.

Reviewed by consultants Per Rydström and Christina Bojesson.

Approved by Oskar Ekelund, consultant and Medical Director at Clinical Microbiology for Kronoberg and Blekinge County.